

# ISO Class 3-5

Aseptic and Sterile Environments

Aseptic and sterile environments isolate the work area from the background environment.

ITW Texwipe has the tools, products and services to maintain the cleanliness of the work area.

## Mini Environment Cleaning and Maintenance Starter Kit

### TX7101 — Mini AlphaMop™ Isolator Cleaning Tool (ICT)

- Economical
- Easier to use and more ergonomic than wipers
  - Cleans hard-to-reach areas quickly and easily
  - Low, flat profile
  - Light weight – can be used with one hand
  - Swivel joint – flexible use
- Excellent for isolators, biosafety cabinets, glove boxes and laminar flow hoods
- Fits into most pass-throughs and ante-chambers
  - Easy to disassemble and reassemble
  - Two handles (18" and 24")
- Tolerates sterilization by autoclave, gamma irradiation or chemical means

### TX3285 — Sterile AlphaSat® 10 with 70% IPA

- Sterile, individually folded, small footprint (13" x 7"), saving valuable space
- Sealed-border wipers prevent fiber release
- Use to wipe down all surfaces, equipment and gloves
- Excellent for spill control in isolators, biosafety cabinets, glove boxes and laminar flow hoods
- Bag-within-a-Bag® cleanroom packaging
  - Double bagging
  - Wipe the exterior bag with isopropyl alcohol (IPA) before taking the inner bag into a cleanroom
- Polyester wipers pre-wetted with 70% IPA
  - Safer and easier to use than a dry wiper and a bottle of IPA
  - Eliminates handling of bottles of flammable liquids
  - Optimal concentration of IPA for killing microbes
- Resealable flexpack
- Use to remove cleaning and disinfectant residues



### TX3211 — SterileWipe™ LP Dry Wiper

- Sterilized for sterile and aseptic processes
- Ideal for cleaning, spill control and wiping gloves
- Bag-within-a-Bag® cleanroom packaging
  - Double bagging
  - Wipe the exterior bag with isopropyl alcohol (IPA) before taking the inner bag into a cleanroom
- Good abrasion resistance
- Tolerates most solvents

### This kit also includes the following literature:

- TechNote ICT-1: Cleaning of Compounding Aseptic Isolators (CAIs) and Class II Biological Safety Cabinets (BSCs)
- Brochures:
  - Products for USP <797> Compliance
  - CleanTips® Swabs
- MSDS: TX3285

For additional technical information go to [www.texwipe.com](http://www.texwipe.com)

	North America	Europe	Asia
Tel	201 684 1800	+45 87 400 220	+65 6468 9433
Fax	201 684 1801	+45 87 400 222	+65 6468 6772
E-mail	<a href="mailto:info@texwipe.com">info@texwipe.com</a>	<a href="mailto:europa@texwipe.com">europa@texwipe.com</a>	<a href="mailto:asia@texwipe.com">asia@texwipe.com</a>

**ITW Texwipe®**

Quality. Consistency. Support.

[www.texwipe.com](http://www.texwipe.com)